

Material description: AMMONIA SOLUTION 35% w/w pure  
CAS No: 1336-21-6

<i>Supplied as:</i>	Product code	Description
	CHE1130	AMMONIA SOLUTION 35% w/w pure 2.5L.
	CHE1132	AMMONIA SOLUTION 35% w/w pure 5L.
	CHE1136	AMMONIA SOLUTION 35% w/w pure 1L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	33			
Residue on evaporation	ppm		25		
Residue on ignition	ppm		10		
Iron; as Fe	ppm		1		
Colour	Hazen		20		
Turbidity	NTU		3		

Revision date: 09/09/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0051/MS/19910909