

Material description: AMMONIUM ACETATE A.R.

CAS No: 631-61-8

*Supplied as:* Product code Description  
CHE1138 AMMONIUM ACETATE A.R. 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97	102		
pH; 5% solution	pH	6.5	7.5		
Chloride	%		0.0005		
Nitrate	%		0.001		
Sulphate	%		0.001		
Iron	%		0.0005		
Lead	%		0.0005		

Revision date: 21/02/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1826/MS/19910221