

Material description: AMMONIUM FERROUS SULPHATE 6H₂O A.R.

CAS No: 7783-85-9

Supplied as: Product code Description
CHE1164 AMMONIUM FERROUS SULPHATE 6H₂O A.R. 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.0			
Copper	ppm		20		
Iron	ppm		200		
Manganese	ppm		500		
Zinc	ppm		30		

Revision date: 10/05/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1958/MS/20210510