

Material description: AMYL ACETATE pure (mixture of isomers)

CAS No: Mixture

Supplied as: Product code Description
CHE1218 AMYL ACETATE pure (mixture of isomers) 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Wt. per ml.	g	0.866	0.874		
Acidity; acetic	%		0.02		
Water	% w/w		0.2		

Revision date: 03/06/1993

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1293/MS/19930603