

Material description: 2-BROMO-2-METHYLPROPANE pure
CAS No: 507-19-7

Supplied as: Product code Description
CHE1248 2-BROMO-2-METHYLPROPANE pure 100ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	96			
Boiling range	°C	71	74		
Wt. per ml.	g	1.215	1.225		

Revision date: 21/02/1991 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1678/MS/19910221