

Material description: BARIUM CHLORIDE 2H₂O pure
CAS No: 10326-27-9

<i>Supplied as:</i>	Product code	Description
	CHE1260	BARIUM CHLORIDE 2H ₂ O pure 500g.
	CHE1262	BARIUM CHLORIDE 2H ₂ O pure 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Insoluble matter	%		0.1		

Revision date: 20/09/2004

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0864/MS/20040920