

Material description: BARIUM SULPHATE FOR SOIL TESTING

CAS No: 7727-43-7

<i>Supplied as:</i>	Product code	Description
	CHE1272A	BARIUM SULPHATE FOR SOIL TESTING 3Kg.
	CHE1274A	BARIUM SULPHATE FOR SOIL TESTING 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Acid insoluble matter	%		0.3		
Soluble Ba compounds	%		0.001		
Acidity/alkalinity; Litmus					Neutral

Revision date: 18/09/1992

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0868/MS/19920918