



Material description: BENZOPHENONE pure
CAS No: 119-61-9

Supplied as: Product code Description
CHE1290 BENZOPHENONE pure 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Melting range	°C	47	49		

Revision date: 01/01/1989 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0094/MS/19890101