

Material description: BICARBONATE INDICATOR SOLUTION
CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE1298	BICARBONATE INDICATOR SOLUTION 100ml.
	CHE1300	BICARBONATE INDICATOR SOLUTION 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Absorbance @ 436nm	AU	0.900	2.000		
Functional test <i>by Colour change</i>					Red to Yellow

Revision date: 21/05/2012 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1891/MS/20120521