

Material description: BIURET REAGENT

CAS No: Mixture

*Supplied as:* Product code

CHE1306

Description

BIURET REAGENT 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
pH @ 20°C	pH	12.8	13.2	13.0	
Density @20°C	g/ml	1.032	1.036	1.034	
Functional test <i>by Colour change</i>					Blue to purple

Revision date: 23/04/2012

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3064/MS/20120423