

Material description: BORIC ACID SOLUTION 4% w/v (with indicator)

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE1314	BORIC ACID SOLUTION 4% w/v (with indicator) 2.5L.
	CHE1316	BORIC ACID SOLUTION 4% w/v (with indicator) 5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
pH 4.0	pH				red-pink
pH 5.0	pH				purple
pH 6.0	pH				green

Revision date: 31/05/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2058/MS/20200531