

Material description: BRIJ 35
CAS No: 9002-92-0

Supplied as: Product code Description
CHE1320 BRIJ 35 100g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Hydroxyl number	mg KOH/g	40	60		
Water	% w/w			3	

Revision date: 13/05/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2288/MS/20200513