



Material description: BROMOCRESOL GREEN SOLUTION

CAS No: Mixture

Supplied as: Product code Description
CHE1336 BROMOCRESOL GREEN SOLUTION 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Yellow to Blue					pH 3.6-5.2

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0895/MS/19890101