



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: BROMOPHENOL BLUE

CAS No: 115-39-9

<i>Supplied as:</i>	Product code	Description
	CHE1340	BROMOPHENOL BLUE 25g.
	CHE1342	BROMOPHENOL BLUE 5g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Yellow to Violet	pH				pH 2.8-4.6

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1897/MS/19890101