



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: BROMOTHYMOL BLUE

CAS No: 76-59-5

Supplied as: Product code Description
CHE1348 BROMOTHYMOL BLUE 5g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Yellow to Blue					pH 6.0-7.6

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1874/MS/19890101