



Material description: BROMOTHYMOL BLUE SOLUTION

CAS No: Mixture

Supplied as:	Product code	Description
	CHE1350	BROMOTHYMOL BLUE SOLUTION 500ml.
	CHE1352	BROMOTHYMOL BLUE SOLUTION 100ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Yellow to Blue	pH				pH 6.0-7.6

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0928/MS/19890101