



Material description: BUFFER TABLETS pH 4.00

CAS No: Mixture

Supplied as: Product code Description
CHE1418 BUFFER TABLETS pH 4.00 100Tabs.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
pH: 1 tablet:100ml DI Water	pH	3.95	4.05		

Revision date: 26/01/2022

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2282/MS/20220126