

Material description: BUTAN-1-OL A.R.

CAS No: 71-36-3

Supplied as: Product code Description
CHE1428 BUTAN-1-OL A.R. 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Boiling range	°C	116	119		
Wt. per ml.	g	0.809	0.813		
Colour	Hazen		10		
Acidity; acetic	%		0.005		
Water	ppm		1000.0		

Revision date: 06/01/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0126/MS/20210106