

Material description: BUTANONE A.R.

CAS No: 78-93-3

<i>Supplied as:</i>	Product code	Description
	CHE1438	BUTANONE A.R. 500ml.
	CHE1440	BUTANONE A.R. 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Initial boiling point	°C			79.1	
Wt. per ml.	g	0.804	0.806		
Colour	Hazen		10		
Acidity; acetic	%		0.003		
Water	% w/w		0.1		
Alcohol impurities	%		0.5		

Revision date: 23/05/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0934/MS/20170523