

Material description: CALCIUM CHLORIDE ANHYDROUS PRILLS

CAS No: 10043-52-4

<i>Supplied as:</i>	Product code	Description
	CHE1484	CALCIUM CHLORIDE ANHYDROUS PRILLS 500g.
	CHE1486	CALCIUM CHLORIDE ANHYDROUS PRILLS 2Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	94			
Sulphate	%		0.1		
Barium	%		0.08		
Iron	ppm		5.00		
Zinc	ppm		2.00		
Copper	ppm		2.00		
Lead	ppm		2.0		
Mercury	ppm		0.030		
Aresnic	ppm		0.030		

Revision date: 08/09/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1530/MS/20210908