

Material description: CHARCOAL WOOD POWDER  
CAS No: 7440-44-0

Supplied as: Product code Description  
CHE1558 CHARCOAL WOOD POWDER 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Particle Size <45um	%	60	85		
Moisture	%		10		
pH	pH	2	4		
Ash	%		4		
Molasses Decolour		500			

Revision date: 19/04/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3334/MS/20180419