

Material description: CHLOROFORM A.R.  
CAS No: 67-66-3

<i>Supplied as:</i>	Product code	Description
	CHE1568	CHLOROFORM A.R. 2.5L.
	CHE1570	CHLOROFORM A.R. 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Water	ppm		500		
Non-volatiles	ppm		10		
1.1-Dichloroethane	ppm		1000		
Bromochloromethane	ppm		2000		
Carbon Tetrachloride	ppm		2000		
Methyl Chloride	ppm		1000		
Chloride	ppm		0.2		
Lead	ppm		0.05		
Organic impurities	ppm		1000		
Free Halogen-C12	ppm		2		
Dichloroethylene	ppm		1000		
Inhibitor-Ethanol	% w/w		0.8		

Revision date: 13/01/2011

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0168/MS/20110113