

Material description: CHROMIUM (III) POTASSIUM SULPHATE 12H<sub>2</sub>O pure  
CAS No: 7788-99-0

*Supplied as:* Product code Description  
CHE1590 CHROMIUM (III) POTASSIUM SULPHATE 12H<sub>2</sub>O pure 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97			
Chromic oxide	%	14.5	15.2		
Iron	%		0.10		

Revision date: 24/09/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0983/MS/19910924