

Material description: CITRIC ACID H2O A.R.

CAS No: 5949-29-1

<i>Supplied as:</i>	Product code	Description
	CHE1600	CITRIC ACID H2O A.R. 500g.
	CHE1602	CITRIC ACID H2O A.R. 2Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5	101.0		
Chloride	%		0.005		
Oxalate	%		0.01		
Lead	%		0.00005		
Arsenic	%		0.0001		
Mercury	%		0.0001		

Revision date: 28/10/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2030/MS/20211028