



Material description: COBALT (II) CHLORIDE 6H<sub>2</sub>O pure  
CAS No: 7791-13-1

*Supplied as:*

Product code	Description
CHE1616	COBALT (II) CHLORIDE 6H <sub>2</sub> O pure 250g.
CHE1618	COBALT (II) CHLORIDE 6H <sub>2</sub> O pure 100g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay (As Co)	% w/w	23.50	25.00		

Revision date: 25/10/2022 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0183/MS/20221025