



**SCIENTIFIC  
LABORATORY  
SUPPLIES**

Material description: CONGO RED (C.I. 22120)

CAS No: 573-58-0

<i>Supplied as:</i>	Product code	Description
	CHE1644	CONGO RED (C.I. 22120) 5g.
	CHE1646	CONGO RED (C.I. 22120) 25g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Violet to Red Orange					pH 3.0-5.0

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2132/MS/19890101