



Material description: COPPER METAL FOIL (0.13mm.)

CAS No: 7440-50-8

*Supplied as:* Product code Description  
CHE1650 COPPER METAL FOIL (0.13mm.) 250g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1700/MS/19890101