

Material description: CRESOL RED
CAS No: 1733-12-6

Supplied as: Product code Description
CHE1662 CRESOL RED 5g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
UV Lambda Max (pH 6.8)	nm	432.0	437.0		
UV Lambda Max (pH 8.8)	nm	571.0	575.0		
Loss on drying @110°C	%		3.0		

Revision date: 20/09/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2975/MS/20210920