



Material description: COPPER (II) CHLORIDE 2H<sub>2</sub>O pure  
CAS No: 10125-13-0

Supplied as: Product code Description  
CHE1692 COPPER (II) CHLORIDE 2H<sub>2</sub>O pure 250g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	96			
Moisture	%		4		

Revision date: 14/05/1992 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0993/MS/19920514