

Material description: COPPER (II) NITRATE 3H<sub>2</sub>O pure  
CAS No: 10031-43-3

<i>Supplied as:</i>	Product code	Description
	CHE1698	COPPER (II) NITRATE 3H <sub>2</sub> O pure 250g.
	CHE1700	COPPER (II) NITRATE 3H <sub>2</sub> O pure 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Chloride	%		0.005		
Sulphate	%		0.02		
Iron	%		0.002		

Revision date: 09/08/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0995/MS/20180809