

Material description: COPPER (II) OXIDE POWDER pure  
CAS No: 1317-38-0

<i>Supplied as:</i>	Product code	Description
	CHE1702	COPPER (II) OXIDE POWDER pure 250g.
	CHE1704	COPPER (II) OXIDE POWDER pure 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97			
Iron	%		0.05		
Acid insoluble matter	%		0.02		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0187/MS/19890101