

Material description: CYCLOHEXANE HPLC  
CAS No: 110-82-7

Supplied as: Product code Description  
CHE1730 CYCLOHEXANE HPLC 2.5L.

## Specification

| Parameter           | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------------|-------|---------------|---------------|---------------|-------|
| Purity              | %     | 99.9          |               |               |       |
| Water               | % w/w |               | 0.02          |               |       |
| UVAbsorption; 210nm | AU    |               | 1             |               |       |
| UVAbsorption; 220nm | AU    |               | 0.3           |               |       |
| UVAbsorption; 250nm | AU    |               | 0.02          |               |       |
| UVAbsorption; 260nm | AU    |               | 0.015         |               |       |
| UVAbsorption; 400nm | AU    |               | 0.01          |               |       |

Revision date: 23/02/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2928/MS/20170223