

Material description: CYCLOHEXANE pure  
CAS No: 110-82-7

<i>Supplied as:</i>	Product code	Description
	CHE1732	CYCLOHEXANE pure 500ml.
	CHE1734	CYCLOHEXANE pure 2.5L.
	CHE1736	CYCLOHEXANE pure 250ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			
Appearance					Clear colourless, free from suspended material
Refractive index @20°C		1.424	1.428		
Water	% w/w		0.05		
Colour	Hazen		10		
Boiling range	°C	80	82		
Wt. per ml.	g	0.774	0.781		

Revision date: 08/06/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1586/MS/20090608