



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: CYCLOHEXANOL pure
CAS No: 108-93-0

<i>Supplied as:</i>	Product code	Description
	CHE1738	CYCLOHEXANOL pure 250ml.
	CHE1740	CYCLOHEXANOL pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	158.5	162.0		
Melting point	°C			22	
Wt. per ml.	g			0.95	

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1005/MS/19890101