

Material description: CYCLOHEXENE pure
CAS No: 110-83-8

Supplied as: Product code Description
CHE1746 CYCLOHEXENE pure 250ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Refractive index @20°C		1.445	1.448		
Infrared Spectrum					Complies

Revision date: 23/01/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1587/MS/20200123