



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: 1,6-DIAMINOHEXANE pure
CAS No: 124-09-4

Supplied as: Product code Description
CHE1776 1,6-DIAMINOHEXANE pure 25g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Melting point	°C	38	44		

Revision date: 28/01/1991 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0214/MS/19910128