

Material description: 4-DIMETHYLAMINOBENZALDEHYDE A.R.

CAS No: 100-10-7

*Supplied as:* Product code Description  
CHE1786 4-DIMETHYLAMINOBENZALDEHYDE A.R. 100g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Melting range	°C	74	75		
Chloride	%		0.005		
Copper	%		0.0005		
Iron	%		0.0005		
Lead	%		0.0005		
Sulphated ash	%		0.05		
Alcohol insolubles	%		0.01		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2393/MS/19890101