

Material description: DICHLOROMETHANE A.R.

CAS No: 75-09-2

<i>Supplied as:</i>	Product code	Description
	CHE1822	DICHLOROMETHANE A.R. 500ml.
	CHE1824	DICHLOROMETHANE A.R. 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.9			
Boiling range	°C	39.5	40.5		
Water	ppm		100		
Non-volatiles	ppm		10		
Acidity; HCl	ppm		5		

Revision date: 22/08/2011

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0238/MS/20110822