

Material description: DICHLOROMETHANE HPLC
CAS No: 75-09-2

Supplied as: Product code Description
CHE1826 DICHLOROMETHANE HPLC 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Acidity/Alkalinity	meq/g		0.00010		
Colour	Hazen		10		
Water	% w/w		0.03		
Non-volatiles	ppm		2.0		
Absorbance @ 233nm	AU		1.000		
Absorbance @ 240nm	AU		0.200		
Absorbance @ 250nm	AU		0.020		
Absorbance @ 260nm	AU		0.005		
Absorbance @ 270nm	AU		0.005		
Absorbance @ 280nm	AU		0.005		

Revision date: 23/02/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2230/MS/20210223