



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: DICHLOROMETHANE pure

CAS No: 75-09-2

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|-----------------------------|
| | CHE1828 | DICHLOROMETHANE pure 2.5L. |
| | CHE1830 | DICHLOROMETHANE pure 500ml. |
| | CHE1832 | DICHLOROMETHANE pure 25L. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Boiling range | °C | 39.0 | 40.5 | | |
| Wt. per ml. | g | 1.323 | 1.327 | | |
| Water | % w/w | | 0.05 | | |
| Non-volatiles | % | | 0.01 | | |

Revision date: 09/05/2018

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0237/MS/20180509