



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: DIETHYL ETHER pure

CAS No: 60-29-7

<i>Supplied as:</i>	Product code	Description
	CHE1842	DIETHYL ETHER pure 2.5L.
	CHE1844	DIETHYL ETHER pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.00			
Boiling range	°C	34	36		
Wt. per ml.	g	0.713	0.718		
Water	% w/w		0.200		

Revision date: 19/02/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0250/MS/20210219