



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: DIETHYLAMINE pure

CAS No: 109-89-7

Supplied as: Product code

CHE1850

Description

DIETHYLAMINE pure 250ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	54	56		
Wt. per ml.	g	0.706	0.708		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1715/MS/19890101