



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: DIMETHYL SULPHOXIDE A.R.

CAS No: 67-68-5

<i>Supplied as:</i>	Product code	Description
	CHE1852	DIMETHYL SULPHOXIDE A.R. 500g.
	CHE1854	DIMETHYL SULPHOXIDE A.R. 2.5Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity-GC	%	99.5			
Non-volatiles	%		0.01		
Titrateable acid, meq/g	%		0.001		
Water	%		0.1		

Revision date: 19/12/2012

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1403/MS/20121219