



Material description: EOSIN Y 24050 (C.I. 45380-water & alcohol soluble)

CAS No: 17372-87-1

Supplied as: Product code Description
CHE1916 EOSIN Y 24050 (C.I. 45380-water & alcohol soluble) 25g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	80	88		
Wavelength of max absorption	nm	515	517		
Specific absorbance		1120	1235		
Loss on Drying	%		8		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3105/MS/19890101