

Material description: ETHANE-1,2-DIOL pure
CAS No: 107-21-1

<i>Supplied as:</i>	Product code	Description
	CHE1918	ETHANE-1,2-DIOL pure 500ml.
	CHE1920	ETHANE-1,2-DIOL pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Assay	%	99.00	101.00		
Boiling range	°C	194	199		
Wt. per ml.	g			1.11	
Acidity; acetic	%		0.02		
Water	% w/w		0.3		

Revision date: 29/04/2021 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1092/MS/20210429