

Material description: ETHANOL ABSOLUTE A.R. (duty paid)

CAS No: 64-17-5

<i>Supplied as:</i>	Product code	Description
	CHE1922	ETHANOL ABSOLUTE A.R. (duty paid) 250ml.
	CHE1924	ETHANOL ABSOLUTE A.R. (duty paid) 2.5L.
	CHE1926	ETHANOL ABSOLUTE A.R. (duty paid) 500ml.
	CHE1928	ETHANOL ABSOLUTE A.R. (duty paid) 100ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.7			
Boiling range	°C	78.0	78.5		
Wt. per ml.	g		0.7897		
Miscibility with water					Complete
Aldehydes	%		0.005		
Non-volatiles	%		0.002		
Water	% w/w		0.2		
Acidity; acetic	%		0.01000		

Revision date: 16/10/2015

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1725/MS/20151016