

Material description: ETHANOL HPLC  
CAS No: 64-17-5

Supplied as: Product code Description  
CHE1930 ETHANOL HPLC 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.7			
Water	% w/w		0.2500		
Absorbance @ 400nm	AU		0.0100		
Absorbance @ 260nm	AU		0.0300		
Absorbance @ 250nm	AU		0.0400		
Absorbance @ 240nm	AU		0.1000		
Absorbance @ 225nm	AU		0.3000		
Absorbance @ 205nm	AU		1.0000		

Revision date: 15/05/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2931/MS/20200515