

Material description: ETHANOL INDUSTRIAL (99% IMS) (74 OP)

CAS No: 64-17-5

<i>Supplied as:</i>	Product code	Description
	CHE1938	ETHANOL INDUSTRIAL (99% IMS) (74 OP) 2.5L.
	CHE1940	ETHANOL INDUSTRIAL (99% IMS) (74 OP) 5L.
	CHE1942	ETHANOL INDUSTRIAL (99% IMS) (74 OP) 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Total alcohols	%	99	100		
Wt. per ml.	g			0.79	

Revision date: 10/08/2015

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0315/MS/20150810