

Material description: ETHANOLAMINE pure
CAS No: 141-43-5

<i>Supplied as:</i>	Product code	Description
	CHE1952	ETHANOLAMINE pure 500ml.
	CHE1954	ETHANOLAMINE pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Water	% w/w		0.5		
Melting range	°C	9.0	11.0		

Revision date: 20/07/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1596/MS/20200720