

Material description: ETHYL ACETATE A.R.

CAS No: 141-78-6

<i>Supplied as:</i>	Product code	Description
	CHE1958	ETHYL ACETATE A.R. 500ml.
	CHE1960	ETHYL ACETATE A.R. 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.7			
Colour	Hazen		10		
Wt. per ml.	g	0.898	0.903		
Water	% w/w		0.10		

Revision date: 02/03/2022

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0318/MS/20220302